

## How to use the Trans-Isolate (T-I) system for isolation and transport of meningococci and other agents causing bacterial meningitis from CSF.

### 1. Procedure for inoculating T-I medium for transporting meningococci and other agents causing bacterial meningitis from CSF:

- 1.1 Remove a vial of Trans-Isolate (T-I) medium from refrigerator at least 30 minutes before inoculating it with the specimen. Allow the vial to warm to room temperature which is more favorable for growth of the organism.
- 1.2 Before inoculating the vial, check to see if there is any visible growth or turbidity. If there is visible growth or turbidity, discard the vial, because it may be contaminated.
- 1.3 Lift up the small lid in the middle of the metal cap on top of the T-I vial.
- 1.4 Disinfect the top of the T-I vial with 70% alcohol or iodine. Allow to dry (usually 30 to 60 seconds).
- 1.5 Use a sterile syringe and sterile needle preferably 21G, 0.8 mm. to aspirate 500 microliters (one-half of an mL) of cerebrospinal fluid (CSF) from the tube containing CSF.
- 1.6 Inject the CSF into the T-I vial through the disinfected dry stopper on the top of the T-I vial.

### 2. Transport and incubation of T-I vials, and inoculation of the culture media

The procedures to follow depend upon how promptly the TI vials can reach the laboratory of reference that will perform culture and isolation.

If T-I vials **cannot** reach the laboratory of reference within 24 hours::

- Label the T-I vial with the date, name of the patient, and any other necessary identifiers.
- Ventilate the T-I vial with a sterile cotton plugged needle. **The** Needle should not dip into the culture media(broth).
- Store the ventilated T-I vial in an upright position at room temperature. Make sure it is away from excessive heat, direct sunlight, and dust.
- Before transporting the vial, remove the ventilating needle from the top of the T-I vial. This will prevent leakage and contamination during shipment.
- Transport the T-I vial in a sealed plastic bag to minimize the risks of contamination and attach the case report form

If TI vials can reach the laboratory of reference within 24 hours:

- Label the T-I vial with the date, name of the patient, and any other necessary identifiers.
- Ship the T-I vials without ventilation.
- Transport the TI in a sealed plastic bag to minimize the risk of contamination and attach the case report form.

**3. Additional recommendations about the proper use of T-I vials and ventilating the inoculated T-I vials:**

- The T-I vials can be used for at least 1 year after the date of production provided that they are stored in the refrigerator.
- Freezing T-I vials destroys the T-I medium.
- Non-inoculated T-I vials should be packed in cold packs for shipment to the laboratory of reference.
- In previous studies (*Ajello et al* below), cultures on ventilated T-I vials 2 to 4 weeks after inoculation with CSF (from patients with acute bacterial meningitis), incubation and transport resulted in a loss of growth in only 20 to 25% of inoculated vials. Without ventilation the losses were much greater.
- Contamination is the single most problematic point with the system. Aseptic measures and understanding the risks are necessary to achieve good recovery of the isolates.

Reference: Ajello GA, Feely JC, Hayes PS, Reingold AL, Bolan G, Broome CV, and Phillips CJ. *Trans-Isolate Medium: a New Medium for Primary Culturing and Transport of Neisseria meningitidis, Streptococcus pneumoniae, and Haemophilus influenzae*. J Clin Microbiol, 1984: 20 (1): 55-58.